

Characterization of Carcino Embryonic Antigen (CEA) Microheterogeneity: A Case Study Considering N-linked Glycosylation Forms



Ravi Amunugama, Michael Ford, Richard Jones and David Allen

MS Bioworks, LLC. 3950 Varsity Drive Ann Arbor, MI 48108




1. Abstract and introduction

Colorectal cancer (CRC) is the fourth leading cause of cancer related deaths in the United States¹. A critical component of CRC care is post-surgical monitoring for cancer recurrence. CEA is a tumor marker for the clinical management of CRC which has the specific utility of monitoring post-operative disease recurrence. After surgery to remove cancerous tissue, the level of CEA in blood can be periodically monitored using an immunoassay. If the levels begin to rise above 6.0 ng/mL there is a high correlation with recurrence of the cancer².

CEA, like many tumor markers is a glycoprotein and there is a significant body of work showing protein glycosylation is greatly affected by diseases such as cancer. We are investigating glycosylation forms of CEA as sensitive and specific biomarkers of CRC. This presentation outlines our methods and includes qualitative and quantitative CEA glycosylation form data.

1. *Cancer Trends Progress Report – 2009/2010 Update*, National Cancer Institute, NIH, DHHS, Bethesda, MD, April 2010, <http://progressreport.cancer.gov>.
2. M.G. Fakih and A. Padmanabhan "CEA Monitoring in Colorectal Cancer" *Oncology* 20, 2006, 1- 15

2. Technology Overview



- 1µg of purified protein (with and without treating with PNGaseF) was separated on a 10% Bis-Tris NuPage gel (Invitrogen) in the MES buffer system.
- The segments corresponding to the modified and unmodified proteins were excised and processed by in-gel digestion using a robot (ProGest, DigiLab) and sequencing grade trypsin (Promega).
- Each digest is processed by LC-MS/MS using a reverse phase gradient.

3. Glycosylated Sites Mapping


IPI00880101: Gene-Symbol = CEACAM5 Protein Amino Acids Sequence (IPI Human Database)

1 XSAPPHRWCI PWQRLLTAS LLTFWNPPT AKLTIESTPF NVAEGKEVLL	N101, N112, N149
51 LVHNLPHLF GYSWYKGERV DGNRQIIGYV IGTQATPGP AYSGREIIVP	N179, N194
101 NASLLIQNI QNDTGFYTLH VIKSDLVNEE ATGQFRVYPE LKPSISINN	N201, N205, N243
151 SKFVLDKDAV AFTCEPETQD ATYLWVWNNQ SLPVSPRLQL SNGNRLTLFL	N253, N271, N285, N289
201 NVTRNDTASY KCETQNPVSA RRSDSVILNV LYGPDAPTIS PLNTSYRSGE	N306, N326, N347
251 NNLNSCHAAS NPPAQYSWV NGTFQSTQE LFIPINIVNN SGSYTCQAIN	N356, N371, N428
301 SDTGLNRTTV TITVYEPK PFITSNNSNP VEDEDAVALT CEPBIQNTTY	N462, N476
351 LWVWVNSQL VSPRLQLSND NRLLTLLSVT RNDVGPYECG IQNKLSDVHS	N504, N525, N549
401 DPVILNVLYG PDDPTISPSY TYRPGVNLV LSCHAASNPP AQYSLWIDGN	N556, N576, N608, N646
451 IQHTQELFI SNITEKNSGL YTCQANNSAS GHSRTTVKTI TVSAELPKPS	N661
501 ISSNNSKPFV DRDAVAFCE PEAQNTTYLW VVNGQSLPVS PRLQLSNGNR	
551 TLTLFNVTRN DARAYVCGIQ NSVSANRSDP VILDVLYGPD TPIISPDPSS	
601 YLSGANLNLV CHSASNSPQ YSWRINGIPQ QHTQVLFIAK ITPNNCTYA	
651 CFVSNLATGR NSIVKSITV SASGTSPLGS AGATVIGIMIG VLVGVVALI	

- 28 potential glycosylation sites are present in CEACAM5 protein (based on NXS/T Sequon).
- Deglycosylation by treating with PNGase F results in diagnostic deamidation of Asparagine (N) residue to Aspartic (D) acid.
- Deamidation of Asn (N) can be readily configured in a database search engine and its presence combined with the motif NXT or NXS can be used as a tool for determining N-linked glycosylation sites.
- A trypsin alone digest is performed to identify background deamidation.
- By this method, nine and/or eleven sites (N194/N549, N201, N205, N371, N476, N201/N556, N576, N646, and N661) were found to be N-glycosylated.


4. Glycan Composition Determination (N194/N549)

Annotated MS2 and MS3 [MS2 of $m/z = 1104.7 (+1)$] Spectra of $m/z = 1225.52 (+2)$ Th. containing glycosylated tryptic peptide LQLSNGNR, Peptide1.



- MS² spectra of Glyco-peptides contain a wealth of information regarding the composition of glycan modification.
- Data are interpreted manually and the m/z ratio of the bare or minimally modified peptide is determined.
- MS³ is performed on the minimally modified peptide to validate primary AA sequence and confirm the N-linked glycosylation site.


5. Heterogeneity of Glycosylation at N194/N549



- Certain other extended glycosylation forms attached to Peptide1 (N194 and/or N549): The figure on left shows glycans without NeuNAc, The figure on right shows glycans with one or more Fucose (Fuc). The glycans containing NeuNAc but no Fucose were also observed (data not shown).
- The carbohydrate compositions indicated above were manually validated by inspection of MS² data.

6. HRAM Quantitation

Extracted Ion Chromatogram (XIC) of m/z 1225.5211



- Peak area values for different glycosylation forms of the same peptide (eg. Peptide1) can be employed for relative quantitation purposes as glycan contribution to the MS response is very minimal.
- High Resolution Accurate Mass data enables sensitive and specific quantitation of individual glycosylation forms in complex backgrounds.


Reprints available from info@msbioworks.com

7. Glycosylation forms identified & validated: Peptide1

Glycosylation form Identifier	Experimental m/z (z)	Composition of the Glycans for Peptide LQLSNGNR	Theoretical m/z (z)	Error (ppm)
1	978.4319 (+2)	(GlcNAc) ₂ (Man) ₃ (Hex)	978.4308 (+2)	1.09
2	998.9454 (+2)	(GlcNAc) ₂ (Man) ₃ (HexNAc)	998.9441 (+2)	1.27
3	1059.4591 (+2)	(GlcNAc) ₂ (Man) ₃ (Hex) ₂	1059.4572 (+2)	1.76
4	1078.4303 (+2)	(GlcNAc) ₂ (Man) ₃ (Hex)(HexNAc)	1079.9705 (+2)	0.71
5	1160.9979 (+2)	(GlcNAc) ₂ (Man) ₃ (Hex) ₂ (HexNAc)	1160.9969 (+2)	0.83
6	1181.5108 (+2)	(GlcNAc) ₂ (Man) ₃ (Hex)(HexNAc) ₂	1181.5102 (+2)	0.48
7	1225.5211 (+2) & 817.3479(+3)	(GlcNAc) ₂ (Man) ₃ (Hex)(HexNAc)(NeuNAc)	1225.5182 (+2)	2.34
8	1242.0250 (+2)	(GlcNAc) ₂ (Man) ₃ (Hex)(HexNAc) ₂	1242.0233 (+2)	1.34
9	1262.5378 (+2)	(GlcNAc) ₂ (Man) ₃ (Hex) ₂ (HexNAc) ₂	1262.5366 (+2)	0.93
10	1306.5474 (+2) & 871.3656 (+3)	(GlcNAc) ₂ (Man) ₃ (Hex) ₂ (HexNAc)(NeuNAc)	1306.5446 (+2)	2.12
11	1327.0606 (+2) & 885.0410 (+3)	(GlcNAc) ₂ (Man) ₃ (Hex)(HexNAc) ₂ (NeuNAc)	1327.0579 (+2)	2.01
12	1387.5736 (+2) & 925.3848 (+3)	(GlcNAc) ₂ (Man) ₃ (Hex) ₂ (HexNAc)(NeuNAc)	1387.5710 (+2)	1.85
13	1408.0877 (+2) & 939.0605 (+3)	(GlcNAc) ₂ (Man) ₃ (Hex) ₂ (HexNAc) ₂ (NeuNAc)	1408.0843 (+2)	2.39
14	952.7355 (+3)	(GlcNAc) ₂ (Man) ₃ (Hex) ₂ (HexNAc) ₂ (NeuNAc)	952.7342 (+3)	1.37
15	987.7462 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	987.7446 (+3)	1.59
16	1006.7537 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1006.7518 (+3)	1.89
17	1020.4288 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1020.4273 (+3)	1.44
18	1055.4401 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc)(NeuNAc)	1055.4378 (+3)	2.21
19	1060.7713 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1060.7694 (+3)	1.79
20	1074.4469 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1074.4449 (+3)	1.83
21	1109.4575 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc)(NeuNAc)	1109.4554 (+3)	1.93
22	1128.4650 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1128.4625 (+3)	2.19
23	1158.1434 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc) ₂ (NeuNAc)	1158.1413 (+3)	1.79
24	1177.1517 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc)(NeuNAc)	1177.1485 (+3)	2.72
25	1182.4828 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1182.4801 (+3)	2.26
26	1225.8369 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc)(NeuNAc)	1225.8345 (+3)	1.99
27	1231.1685 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc)(NeuNAc)	1231.1661 (+3)	1.95
28	1250.1761 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1250.1733 (+3)	1.47
29	1279.8515 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc) ₂ (NeuNAc)	1279.8521 (+3)	0.44
30	1298.8609 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc)(NeuNAc)	1298.8582 (+3)	1.29
31	1304.1937 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (NeuNAc)	1304.1909 (+3)	2.18
32	1328.5360 (+3)	(GlcNAc) ₂ (Man) ₃ (HexNAc) ₂ (Hex) ₂ (Fuc) ₂ (NeuNAc)	1328.5380 (+3)	1.53

- Fifty (50) Multiple glycoforms of the LQLSNGNR peptide are identified from a single purified protein sample.
- High resolution accurate mass data are essential to support the identification and characterization process

8. Glycosylation Forms Quantification



- These data show the relative abundance of the glycosylation forms (1-32) identified in Section 7.0.
- Based on these data the (GlcNAc)₂(Man)₃(Hex)₂ glycosylation form of LQLSNGNR is the most abundant.
- Under close analysis, there is clear variation in relative abundance data that is a function of the carbohydrate composition of the glycans.

9. Summary

- This presentation summarizes a workflow for the identification, characterization, and quantification of protein glycosylation forms.
- The combination of immuno-purification, a targeted multiple enzyme approach and high resolution accurate mass data enables the accurate identification of N-glycosylation sites and their glycosylation forms
- Glycosylation screening with PNGaseF identified nine sites as N-linked glycosylated.
- Greater degree of heterogeneity of glycosylation observed for N194/N549 (50 forms) and N201/N556 (43 forms), out of five/seven sites characterized so far. Other sites include N476 (12 forms), N371 (17 forms) and N576 (13 forms).
- The glycans which are sialylated, contain only one NeuNAc (at least in those validated so far).
- These data represent only a small portion of the quantitative data that has been derived for the different glycosylation forms of CEA and have value not only for understanding glycan biosynthesis but also for understanding the value of protein micro-heterogeneity in evaluating the disease state of CRC patients.